




UTI COMBO STRIPS



UTI Combo Strips test for white blood cells (Leukocytes) and Nitrite in urine.

The presence of specific compounds in urine is used to determine the approximate values of the various analytes. Reactions between chemicals on the pads and said compounds in urine, combined with certain color indicators, produce colors that correlate with the presence and concentration of the relevant analyte.



UTI Combo Strips
REF #52301 for 1 test

LOOK AT THESE GREAT FEATURES!



Provides quick & easy results in 2 minutes.



One test. Two results!



Same test procedure as the leading national brand!



510(k) cleared for OTC use.



Store at room temperature.

Detect Urinary Tract Infections In 2 Minutes!



Leukocytes

Detects Leukocytes levels as low as 15 WBC's/ μ L.



Nitrite

Detects Nitrite levels as low as 0.075 to 0.10mg/dL.



www.germainelabs.com



800.854.8446
210.692.4192
210.692.4198 Fax

THE TESTS YOU WANT.
THE ACCURACY YOU NEED!

DISTRIBUTED BY:

